



**Be Right™**



## **Funnel, analytical, 65 mL approx. volume**

<b>Product #:</b>	<b>108367</b>
AED Price:	Contact Hach
Available	

Standard long stem funnel meets the needs of a wide variety of applications. Mounded ribbing prevents "air lock" commonly occurring during funneling and filtrations. A steep, 58-degree cone slope makes for fast and easy flow of either liquids or small solids. Polypropylene construction is resistant to chemicals and makes this funnel fully autoclavable for sterilization.

---

### **Specifications**

Quantity:	1 each
-----------	--------