



Tube, Culture, 16 x 100 mm, 396/pk

| | |
|-------------------|----------------|
| Product #: | 2899300 |
| AED Price: | Contact Hach |
| Available | |

Caps (**SOLD SEPARATELY** PN 2241106) have a firmly bonded silicone rubber liner, providing a sure and highly inert sealing surface.

Specifications

| | |
|-----------|------------------------|
| Capacity: | Capacity to Neck 11 mL |
| Diameter: | 16 mm |
| Length: | 100 mm |
| Quantity: | 396 pk |