



## Bottle, sampling, polycarbonate, 125 mL, 6/pk

|            |              |
|------------|--------------|
| Product #: | 2324333      |
| AED Price: | Contact Hach |
| Available  |              |

Transparent polycarbonate bottle with cap features wide mouth (3.8cm/1.48" diameter) opening for easy filling and square design for maximum space efficiency. Autoclavable with molded-in graduations (do not engage threads for autoclave cycle.)

---

### Specifications

|                   |         |
|-------------------|---------|
| Body Material:    | Plastic |
| Capacity:         | 125 mL  |
| Material:         | PC      |
| Package size UOM: | 6pk     |
| Quantity:         | 6 /pk   |