



Cylinder, graduated, mixing, glass, 25 mL ± 0.3 mL, 0.5 mL divisions, polyethylene stopper

Product #:	2088640
AED Price:	Contact Hach
Available	

This glass graduated cylinder is calibrated “to contain” with a molded in base. Each includes a polyethylene stopper. Graduations are in durable blue enamel bake.

Specifications

Capacity:	25 mL
Certifications:	Class B
Housing material:	HDPE
Measurement method:	To Contain
Quantity:	1 each
Type:	Single Metric, Mixing