



Buret, Class A, Glass, 25 mL

| | |
|----------------------|--------------|
| Product #: | 2636540 |
| AED Price: | Contact Hach |
| Ships within 2 weeks | |

Glass, funnel-fill buret. Certified Class A (ASTM Specification E287) with permanent serial number and certificate for NIST traceability. PTFE stopcock does not require grease. Dust cap included. 10 mL buret meets requirements for potentiometric titrations (ASTM D664).

Specifications

| | |
|-------------------|---------------|
| Capacity: | 25 mL |
| Housing material: | Class A Glass |
| Special Features: | Buret |