



Cylinder, graduated, 1000 mL, w/ handles

Product #:	2612953
AED Price:	Contact Hach
Available	

One-piece polypropylene construction with a molded in handle top ensure a secure, solid grip and superior control for safer, more convenient pouring. The “Total Contained, Total Delivered ”double scale is accurate and easy to read (cylinder will not create a meniscus), while the hexagonal base prevents tipping or rolling. Hand opening is 89x 29 mm (3 x 11/8"). Steam autoclavable upto 121°C (250°F).

Specifications

Capacity:	1000 mL
Certifications:	None
Housing material:	PP
Measurement method:	To Deliver/To Contain
Quantity:	1 each
Type:	Double Metric