



Cylinder, graduated, mixing, Glass, 50 mL \pm 0.4 mL, 1.0 mL divisions, glass stopper #16

Product #:	2088641
AED Price:	Contact Hach
Available	

This tall form glass graduated cylinder is calibrated “to contain” with a molded in base. Each includes glass penny head stopper. Graduations are in durable white enamel bake.

Specifications

Capacity:	50 mL
Certifications:	Class B
Housing material:	Glass
Measurement method:	To Contain
Quantity:	1 each
Type:	Single Metric, Mixing