



## Cylinder, Graduated, 500 mL, Certified

<b>Product #:</b>	<b>2636249</b>
AED Price:	Contact Hach
Ships within 2 weeks	

Certified Class A, this glass cylinder is calibrated "to deliver" according to ASTM E542 & E1272. Serialized and traceable to NIST standards (includes Certificate of Identification and Capacity). Single metric scale in white marking. Bumper guard included for all sizes  $\geq 25$  mL.

---

### Specifications

Capacity:	500 mL
Certifications:	Class A
Height:	390 mm
Housing material:	Glass
Measurement method:	To Deliver
Quantity:	1 each
Type:	Single Metric