



## Cylinder, graduated, 100 mL +- 0.6 mL, 1.0 mL divisions (white markings)

<b>Product #:</b>	<b>50842</b>
AED Price:	Contact Hach
Available	

Calibrated "to deliver" glass cylinder. Single metric scale in durable white markings. Bumper guards included for all sizes  $\geq$  25 mL.

---

### Specifications

Capacity:	100 mL
Certifications:	Class B
Height:	255 mm
Housing material:	Glass
Measurement method:	To Deliver
Quantity:	1 each
Type:	Single Metric