



Cylinder, Graduated, 250 mL, Certified

Product #:	2636246
AED Price:	Contact Hach
Ships within 2 weeks	

Certified Class A, this glass cylinder is calibrated "to deliver" according to ASTM E542 & E1272. Serialized and traceable to NIST standards (includes Certificate of Identification and Capacity). Single metric scale in white marking. Bumper guard included for all sizes \geq 25 mL.

Specifications

Capacity:	250 mL
Certifications:	Class A
Height:	340 mm
Housing material:	Glass
Measurement method:	To Deliver
Quantity:	1 each
Type:	Single Metric