



Sample Container, 120 mL High PF

Product #:	2552200
AED Price:	Contact Hach
Available	

These 120 mL (4oz.) High Profile Corning Snap-Seal™ polypropylene containers are designed for a wide variety of applications and provide a reliable leak-tight seal when closed properly. Approximately 3.5 inches tall (91 mm) and 1.75 inches (45 mm) diameter. Plastic, snap seal, translucent. Sold individually.

Specifications

Capacity:	120 mL
Diameter:	45 mm (approx)
Quantity:	1 /pk