



## Cylinder, graduated, polypropylene, 100 mL $\pm 0.6$ mL, 1.0 mL divisions

|                   |               |
|-------------------|---------------|
| <b>Product #:</b> | <b>108142</b> |
| AED Price:        | Contact Hach  |
| Available         |               |

---

### Specifications

|                     |                       |
|---------------------|-----------------------|
| Capacity:           | 100 mL                |
| Certifications:     | None                  |
| Housing material:   | PP                    |
| Measurement method: | To Deliver/To Contain |
| Quantity:           | 1 each                |
| Type:               | Single Metric         |