



## Bottle stoppers BOD, acrylic, 25/pk

Product #:	2943900
AED Price:	Contact Hach
Available	

### Disposable BOD Bottle

Plastic stoppers for 60 and 300 mL glass or plastic BOD bottles. 25/pk.

**Fresh bottle, ready for use**

**No cleaning required**

**Unbreakable and recyclable**

**Oxygen will not diffuse in or out of sample**

**EPA compliant**

---

### Specifications

Capacity:	300 mL
Method number:	8043
Quantity:	25 /pk
Units:	25/pk