



Parafilm M Laboratory Film

Product #:	251764
AED Price:	Contact Hach
Ships within 2 weeks	

Seals labware quickly, effectively. Film prevents sample evaporation and contamination; protects anhydrous materials. Can be folded repeatedly and twisted around sharp surfaces without rupturing. Stretched lengthwise at 21 °C (71 °F), Parafilm M increases approximately 200%. Becomes more flexible, softer, stronger as temperature increases. Resistant to alcohol, medium concentration acid, and alkalies. Width: 10 cm (4"); 38 m (125') roll.

Specifications

Quantity:	1 each
-----------	--------