



## Test tube, rack, polyethylene, 16 mm (O.D.), 60 holes

|                   |                |
|-------------------|----------------|
| <b>Product #:</b> | <b>2497901</b> |
| AED Price:        | Contact Hach   |
| Available         |                |

Sixty-place rack features a special blend of polypropylene which sinks in water, assuring non-floating racks when immersed. Autoclavable. Stackable and autoclavable. Overall dimensions: 24.6 x 10.7 x 6.4 cm (9.7 x 4.2 x 2.5"). Color may vary.

---

### Specifications

|                        |        |
|------------------------|--------|
| Diameter:              | 16 mm  |
| Quantity:              | 1 each |
| Test tube arrangement: | 5 x 12 |
| Test tube holes:       | 60     |