



## Stir bar, micro, 2 x 7 mm

|                   |                |
|-------------------|----------------|
| <b>Product #:</b> | <b>2054960</b> |
| AED Price:        | Contact Hach   |
| Available         |                |

Use micro stir bar for mixing in vials, test tubes, spectrophotometer cells and other small-volume containers. Usable at temperatures up to 260°C. Alnico V magnet encapsulated in PTFE.

---

### Specifications

|           |        |
|-----------|--------|
| Quantity: | 1 each |
|-----------|--------|