



Test kit, boiler feed & scale, model BSC-1, drop count titration

Product #:

2350300

AED Price:

Contact Hach

Ships within 2 weeks

Hazardous



Proper maintenance and economy of operation of modern boilers and cooling towers require a careful balance between water quality parameters and treatment chemicals. An economical test kit for monitoring scale formation in small industrial boilers, the BSC-1 Boiler Feed and Scale Test Kit uses drop count titration methods to determine levels of alkalinity, chloride, and hardness.

Specifications

Application:	Boiler Water
Model:	BSC-1
Platform:	Drop Count Titrator