



Iron and manganese color disc test kit, model IR-20

Product #:

146300

AED Price:
Available

Contact Hach

Hazardous



When customers complain of stained and dingy laundry, this kit lets you do quick, on-the-spot tests for iron and manganese to see where the problem lies.

Method / Chemistry: Color disc colorimetry. Test: Iron/ Manganese. Range (mg/L): 0-5/ 0-3.

Simple, fast, color disc method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifications

| | |
|----------------------------|--------------------|
| Application: | Water Quality |
| Measuring range: | 0 - 3.0 |
| Measuring range (2): | 0 - 4.0 mg/L Fe |
| Model: | IR-20 |
| Parameter: | Iron and Manganese |
| Platform: | Color Disc |
| Smallest Increments Steps: | 0.2 |