



Test kit, pool master

Product #:

2315001

AED Price:
Available

Contact Hach

Hazardous



Hach's Pool Master Test Kit combines several tests into a single test kit. The test kit includes tests for alkalinity, chlorine, hardness, and pH, plus specialized tests for total alkalinity, calcium hardness, and cyanuric acid.

Uses simple, economical drop count titration and color disc methods

Specifications

Application:	Pool and Spa
Measuring range:	6.6 - 8.4 pH units
Method:	pH
Model:	Pool Master
Platform:	Color Disc