



## Iron Test Kit, Low Range (Model IR-24)

**Product #:**

**255600**

**Hazardous**

**AED Price:**

**Contact Hach**

**Ships within 2 weeks**



Simple, low-cost colour disc method for Iron determination. The Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field.

---

### Specifications

Case style:	D
Measuring range:	Dual Range
Measuring range (2):	0 - 1.0 mg/L Fe
Method:	Colour disc/FerroZine
Model:	IR-24
Number of tests:	100
Parameter:	Iron Ferrous
Platform:	Color Disc
Smallest increment:	0.002
	0.01
Standard range:	0 - 0.2 mg/L Fe; 0 - 1.0 mg/L Fe
Test kit model:	IR-24
Type:	Double Metric