



Test kit, nitrogen, ammonia, model NI-SA

Product #:

2428700

Hazardous

AED Price:

Contact Hach

Ships within 3 days



Designed for testing ammonia in freshwater or seawater, this kit's color disc method offers low-cost, fast, reliable results in an easy-to-carry kit.

Designed for testing ammonia in freshwater or seawater

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Specifications

Case style:	D
Measuring range:	0 - 2.4 mg/L NH ₃ -N
Method:	Color disc/Salicylate
Number of tests:	100
Parameter:	Nitrogen, Ammonia, mid range - For freshwater or seawater**
Platform:	Color Disc
Smallest increment:	0.1
Test kit model:	NI - SA