



**Be Right™**



## **Bottle, Sampling, Polycarbonate, 500 mL, 24/pk**

<b>Product #:</b>	<b>2324579</b>
AED Price:	Contact Hach
Available	

Transparent polycarbonate bottle with cap features wide mouth (3.8 cm/1.48" diameter) opening for easy filling and square design for maximum space efficiency. Autoclavable with molded-in graduations (do not engage threads for autoclave cycle.)

---

### **Specifications**

Body Material:	Polycarbonate
Capacity:	500 mL
Quantity:	24 /pk