



Nitrate standard, 200 mg/L NO₃ (45.2 mg/L NO₃-N), 500 mL

Product #:	LCW827
AED Price:	Contact Hach
Available	

Standard solution for online nitrate determination

Specifications

Concentration:	200 mg/L NO ₃ (45.2 mg/L NO ₃ -N)
Instrument:	Nitratax plus sc, Nitratax clear sc, NT3100sc, NT3200sc
Parameter:	Nitrate
Platform:	Solution
Quantity:	500 mL

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)