



Be Right™



Copper Standard Solution, 10 mg/L Cu (NIST), 100 mL MDB, for APA6000

Product #:	12932
AED Price:	Contact Hach
Available	

Specifications

Concentration:	10 mg/L
Description:	Copper Standard Solution as Cu (NIST)
Instrument:	APA6000
Parameter:	Hardness
Quantity:	100 mL
Storage conditions:	Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture.
Volume:	100 mL MDB

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)