



pH Buffer Solution, pH 7.00, 500 mL

Product #:	2283549.99
AED Price:	Contact Hach
Available	

Quality Buffer Solution, colour-coded yellow for easy identification. NIST traceable. pH is 7.00 ± 0.02 at 25°C.

Specifications

Color coded?:	Yes
Concentration:	7.00 ± 0.02 pH
Description:	Quality Buffer Solution, colour-coded yellow (NIST)
Package Type:	Bottle
pH value(s):	7.00
Platform:	Solution
Quantity:	500 mL