



## COD TNT, 0 - 15000 mg/L O<sub>2</sub>

**Product #:** 2415951  
**AED Price:** Contact Hach

Hazardous



Reagent Vials (Ø 16 mm) contain premeasured COD reagent, optimised for reliable measuring results and easy handling. TNT = Test 'N Tube

---

### Specifications

|                     |                                    |
|---------------------|------------------------------------|
| Instrument:         | DR900                              |
| Measuring range:    | 200 - 15000 mg/L O <sub>2</sub>    |
| Method:             | Dichromate                         |
| Method number:      | 8000                               |
| Number of tests:    | 25 (1 cuvette for blank needed)    |
| Parameter:          | COD                                |
| Platform:           | Test 'N Tube                       |
| Storage conditions: | 15 °C - 25 °C (protect from light) |

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)