



## Sension+ pH combination puncture electrode for food applications, for portable meters

**Product #:** LZW5053.97.0002

**AED Price:** Contact Hach

pH puncture electrode for portable meters, solid gel filled, encapsulated Ag/AgCl reference, 1 ceramic and 1 open diaphragm. Glass electrode, 316 stainless steel body. Range 2-14 pH, 0-60 °C. 1 m cable, MP5 connector.

The electrode is suitable for food applications, e.g cheese, meat, fish, fruits, bread doughs, etc.

**316 stainless steel body**

**Encapsulated Ag/AgCl reference.**

---

### Specifications

Accuracy:	0.02 pH with pH1 Meter
Cable length:	1 m
Connector:	Sension+ MP5
Electrode junction :	1 ceramic, 1 open diaphragm
Electrode type:	Portable
Filling solution:	Non-refillable gel
Material:	External body: 316 stainless steel, pH sensor: Glass
Measuring range:	2 - 14 pH
Parameter:	pH
Sensor type :	Glass
Special Features:	Puncture probe, stainless steel body.
Temperature range:	0 - 60 °C
Type:	Portable
Warranty:	6 months