



## Thermostat cabinet for BOD Direct Plus

**Product #:** LZQ080  
**AED Price:** Contact Hach

Thermostat-controlled incubator specially designed for BOD determination at 20 °C.

Allows continuous regulation of a temperature range of 2 - 40 °C for various fields of application. The temperature can be set in increments of 0.1 °C.

Equipped with four standard 230 Volt earthed sockets inside, a low power consumption and LED display for temperature.

---

### Specifications

Compliance certifications:	CE
Dimensions (H x W x D):	850 mm x 600 mm x 600 mm
Manual languages:	DE/EN/FR/ES
Power supply:	220 - 240 V / 50 Hz
Warranty:	24 months
Weight:	39 kg
Width:	mm