



**Be Right™**



## High Range Orthophosphate Chemkey Reagents (box of 25)

**Product #:** 8636500

**AED Price:** Contact Hach

**Use with the SL1000 Portable Parallel Analyser (PPA) for dramatically streamlined water quality testing.**

High Range Orthophosphate Chemkeys, range 2.0 - 30.0 mg/L PO<sub>4</sub>. For SL1000 Portable Parallel Analyser (PPA).

**Current SL1000 instrument customers who want to test Orthophosphate need to perform a software update that can be located on the instrument download section.**

---

### Specifications

|                  |                                 |
|------------------|---------------------------------|
| Measuring range: | 2.0 - 30.0 mg/L PO <sub>4</sub> |
| Parameter:       | Orthophosphate                  |
| Platform:        | Chemkeys                        |
| Quantity:        | pk/25                           |