



Low Range Orthophosphate Chemkey Reagents (box of 25)

Product #: 8636600
AED Price: Contact Hach

Use with the SL1000 Portable Parallel Analyser (PPA) for dramatically streamlined water quality testing.

Low Range Orthophosphate Chemkeys, range 0.20 - 4.00 mg/L PO₄. For SL1000 Portable Parallel Analyser (PPA).

Current SL1000 instrument customers who want to test Orthophosphate need to perform a software update that can be located on the instrument download section.

Specifications

Measuring range:	0.20 - 4.00 mg/L PO ₄
Parameter:	Orthophosphate
Platform:	Chemkeys
Quantity:	pk/25