



Be Right™



High Range Total Alkalinity Chemkey Reagents (box of 25)

Product #:

8636100

AED Price:

Contact Hach

Use with the SL1000 Portable Parallel Analyser (PPA) for dramatically streamlined water quality testing.

High Range Total Alkalinity Chemkeys, range 200 - 700 mg/L CaCO₃. For SL1000 Portable Parallel Analyser (PPA).

Current SL1000 instrument customers who want to test Alkalinity need to perform a software update that can be located on the instrument download section.

Specifications

Measuring range:	200 - 700 mg/L CaCO ₃
Method:	Colorimetric
Method number:	10283
Parameter:	Alkalinity, HR
Platform:	Chemkeys
Quantity:	pk/25