



## Sulfuric Acid Standardette Solution #100 pk 12.

**Product #:**

**ATS100**

**Hazardous**

**AED Price:**

Contact Hach



### Save time and get more accurate titration results

The Sulfuric Acid Standardette Solution #100 was designed to ensure precision, efficiency and standardization of 0.100 N base titrants when performing titrations in your lab. With a unique method of standardizing solutions, these single use bottles will make your process faster, easier and more accurate than traditional methods.

These Standardettes are used for titrations across many industries and are great for use with the Hach AT1000 product line.

---

### Specifications

Concentration: 0.100 N

What's included?: Box of 12 Standardettes, each containing Sulfuric Acid for the standardization 0.100 N Base. For use with the AT1000.

---

### What's included?

Box of 12 Standardettes, each containing Sulfuric Acid for the standardization 0.100 N Base. For use with the AT1000.