



Ammonia, Free & Total, and Monochloramine Chemkey Reagents (box of 25 each)

Product #: 8791900
AED Price: Contact Hach

Ammonia, free 0.05 - 0.5 mg/L $\text{NH}_3\text{-N}$ / total 0.05 - 1.50 mg/L $\text{NH}_3\text{-N}$, and Monochloramine, 0.04 - 4.0 mg/L Cl_2 , Chemkeys. For SL1000 Portable Parallel Analyser (PPA).

Specifications

Kit Option:	Free & Total Ammonia, and Monochloramine
Measuring range:	Free Ammonia: 0.05 - 0.50 mg/L $\text{NH}_3\text{-N}$ Total Ammonia: 0.05 - 1.50 mg/L $\text{NH}_3\text{-N}$ Monochloramine: 0.04 - 4.00 mg/L Cl_2
Parameter:	Ammonia, Free & Total; Monochloramine
Platform:	Chemkeys
Quantity:	pk/50