



## **Titration Application Pack Beverage: pH, Total Acidity in Beer, Kombucha, Soft Drink/Juice/Energy Drink Bottlers, Single 25 mL**

**Product #:** AP0025.AT1102

**AED Price:** Contact Hach

Application Pack for the determination of pH, Total Acidity in Beer, Kombucha, Soft Drink/Juice/Energy Drink Bottlers. Suited for Titrallab AT1000 Series Automatic Titrator Models AT1102/AT1112/AT1122/AT1222.

---

### **Specifications**

Application:	Food & Beverage
Compatible instruments:	AT1102
	AT1112
	AT1122
	AT1222
Parameter:	pH, Total Acidity (Single Syringe)
Syringe:	Single
Type:	Application Pack
Warranty:	24 months
What's included?:	Application pack contains a pH probe, USB key with pre-programmed applications, polypropylene beakers, magnetic stirring bars, bottle stoppers, conical adapter and one (1) high precision glass 25 mL syringe with syringe ring.

The AT1000 Instrument and Reagents are sold separately.

---

### **What's included?**

Application pack contains a pH probe, USB key with pre-programmed applications, polypropylene beakers, magnetic stirring bars, bottle stoppers, conical adapter and one (1) high precision glass 25 mL syringe with syringe ring. The AT1000 Instrument and Reagents are sold separately.