



Peracetic Acid Chemkey Reagents (box of 25)

Product #: 8635200
AED Price: Contact Hach

Peracetic Acid Chemkeys, range 0.04 - 50.0 mg/L PAA. For SL1000 Portable Parallel Analyzer (PPA).

Specifications

| | |
|-------------------|----------------------|
| Instrument: | SL1000 |
| Measuring range: | 0.04 - 50.0 mg/L PAA |
| Method number: | 10297 |
| Parameter: | Peracetic Acid |
| Platform: | Chemkeys |
| Quantity: | pk/25 |
| What's included?: | 25 Chemkeys |

What's included?

25 Chemkeys