



pH combination electrode with TOP68 connector, PT100, PG 13.5

Product #: LZX546.1
AED Price: Contact Hach

Sterilisable process electrode and double reference system for higher temperature and pressure, with head connector TOP68.

Specifications

Cable connection:	TOP68
Diameter:	12 mm
Diaphragm:	Ring diaphragm, porous PTFE, double
Electrolyte:	KCl/AgCl + KNO ₃ gel, suitable for high temperatures
Length:	120 mm
Measuring range:	0 - 14 pH
Operating temperature range:	0 - 135 °C
Pressure range:	16 bar
Reference system:	Ag/AgCl
Temperature sensor:	PT100
Warranty:	6 months
Weight:	Approximately 0.15 kg