



## Ammonium cuvette test 100-1,800 mg/L $\text{NH}_4\text{-N}$ , 25 tests

Product #:

LCK502

AED Price:

Contact Hach

Hazardous



### LCK cuvette tests: Minimum amounts of reagents, maximum safety.

Pre-dosed reagents for maximum safety. Easy to use and photometric evaluation, packaged in a fully equipped box. Environmentally friendly due to less chemicals.

Please note that a 20  $\mu\text{L}$  micropipette (2946003) is required for this test.

#### Precise and reliable measurement values

Maximum safety for users, thanks to the closed system and low amounts of reagents

Convenient and error-free dosing of the reagents

Barcode label for automatic recognition in the photometer

Differentiation between tests and measuring ranges by means of colour-coding

---

### Specifications

According to standard:	ISO 7150-1, DIN 38406 E5-1, UNI 11669:2017, ISO 23695
Instrument:	DR3900, DR6000, DR1900
Measuring range:	100 - 1800 mg/L $\text{NH}_4\text{-N}$
Method:	Indophenol blue
Number of tests:	25
Parameter:	Ammonium
Platform:	LCK
Storage conditions:	2 - 8 °C (keep refrigerated)