



## Nitrite cuvette test 2 - 90 mg/L NO<sub>2</sub>-N, 25 tests

Product #:

LCK343

Hazardous

AED Price:

Contact Hach



### LCK cuvette tests: Minimum amounts of reagents, maximum safety.

Pre-dosed reagents for maximum safety. Easy to use and photometric evaluation, packaged in a fully equipped box. Environmentally friendly due to less chemicals.

**Precise and reliable measurement values**

**Maximum safety for users, thanks to the closed system and low amounts of reagents**

**Convenient and error-free dosing of the reagents**

**Barcode label for automatic recognition in the photometer**

**Differentiation between tests and measuring ranges by means of colour-coding**

---

### Specifications

According to standard:	EN ISO 26777, DIN 38405 D10
Instrument:	DR3900, DR6000, DR1900
Measuring range:	2 - 90 mg/L NO <sub>2</sub> -N
Method:	Diazotization
Number of tests:	25
Parameter:	Nitrite
Platform:	LCK
Storage conditions:	15 - 25 °C