



## CL17sc Cell Cleaning Kit

**Product #:**

**8573100**

**AED Price:**

Contact Hach

**Hazardous**



Hach recommends cleaning the CL17sc cell on a monthly (or as needed) basis. Frequency of cleaning may depend on the application and sample qualities.

The CL17sc Cell Cleaning Kit contains 5.25N sulfuric acid (100 mL dropper bottle) and cotton swabs (10).

---

### Specifications

Instrument:	CL17sc
Measuring range:	N/A
Parameter:	N/A
Storage conditions:	Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture.
What's included?:	Sulfuric acid (100 mL dropper bottle) and cotton swabs (10).

---

### What's included?

Sulfuric acid (100 mL dropper bottle) and cotton swabs (10).