



## CL17sc Calibration Verification Kit

<b>Product #:</b>	<b>8568200</b>
AED Price:	Contact Hach

The CL17sc Calibration Verification Kit includes two bottles of deionized water (500 mL), a chlorine standard (20 mL ampule), ampule breaker, syringe (50 mL), tubing, and fittings.

The CL17sc calibration curve is set at the factory for performance to specifications. User adjustment of the factory calibration curve is only recommended in two instances: if required by a regulatory agency for compliance reporting purposes or if a large repair of the analyzer is done.

---

### Specifications

Instrument:	CL17sc
Measuring range:	N/A
Parameter:	N/A
What's included?:	2x 500 mL bottle of deionized water, 20 mL ampule of chlorine standard, ampule breaker, 50 mL syringe, tubing, and fittings.

---

### What's included?

2x 500 mL bottle of deionized water, 20 mL ampule of chlorine standard, ampule breaker, 50 mL syringe, tubing, and fittings.