



Be Right™



Bottle Set of 12, 2 L, Polyethylene

Product #:	9495600
AED Price:	Contact Hach

Bottle Set, 12 - 2 Liter Polyethelyne with Caps

Available for SD900 Refrigerated and All Weather Samplers

Specifications

What's included?:	#12 ea. 2 Liter Polyethylene Bottles p/n 9493000
	#12 ea. caps p/n 621

What's included?

#12 ea. 2 Liter Polyethylene Bottles p/n 9493000#12 ea. caps p/n 621