



Orthophosphate Standard Solution 0.2 mg/L $\text{PO}_4\text{-P}$

Product #:	LCA549
AED Price:	Contact Hach

Orthophosphate Standard Solution with 0.2 mg/L $\text{PO}_4\text{-P}$ for quality control checks with LCK549 Orthophosphate cuvette test. LCA549 is also suitable for quality checks with Phosphax sc LR.

Specifications

Concentration:	0.2 mg/L $\text{PO}_4\text{-P}$
Parameter:	Orthophosphate
Quantity:	100 mL
Storage conditions:	15 - 25 °C