



## Portable HQ Series Standard Field Case for Standard Probes

**Product #:** LEZ015.99.A001A  
**AED Price:** Contact Hach

Rugged, Convenient Storage for Field Use

Rugged Field Case features customer molded case to hold portable meter, a maximum of 3 standard probes with up to 3 meters of cable, and accessories in place.

All-in-One durable rugged field cases provide convenience and include space for sample collection containers, meter and standard probes, along with flexible space for other field accessories.

---

### Specifications

**What's included?:** Includes Standard Field Case, three sample cups, and removable foam insert

---

### What's included?

Includes Standard Field Case, three sample cups, and removable foam insert