



## Fluoride Chemkey Reagents (box of 25)

**Product #:** 9878000  
**AED Price:** Contact Hach

### Fluoride Chemkey for SL1000 Portable Parallel Analyzer

For determination of Fluoride using Hach Method 10312

Ensure SL1000 has been updated with Software Version V3.0.12.1489 to use this Chemkey.

#### Less Variability

Automation and internal temperature control make the entire process consistent and repeatable, while applying the same processes as current Hach methods. Avoid manual steps that can introduce variability, even when performed by experienced testers.

#### Less Hassle

All chemicals and processes are entirely contained inside the Chemkey. There are no powder pillows or glass vials to handle.

---

### Specifications

|                  |                   |
|------------------|-------------------|
| Instrument:      | SL1000            |
| Measuring range: | 0.10 - 4.0 mg/L F |
| Number of tests: | 25                |
| Parameter:       | Fluoride          |