



pH Buffer Solution, pH 10.01, 1000 mL

Product #:	2283653
AED Price:	Contact Hach
Available	

pH Buffer Standard Solution for pH electrode calibration. Colour-coded blue pH Buffer solution for easy identification. pH is 10.01 ± 0.02 at 25°C.

Specifications

Color coded?:	Yes
Concentration:	10.01 pH \pm 0.02 pH at 25 °C
Parameter:	pH
pH value(s):	10.01 pH \pm 0.02 pH at 25 °C
Platform:	Solution
Quantity:	1000 mL