



ULR Free Chlorine Fluorescence Test Kit, 2-100 µg/L (ppb), 100 Tests

Product #:

34251000

Hazardous

AED Price:

Contact Hach



Free Chlorine Test Kit for use with the DR1300 FL Fluorometer.

Specifications

Measuring range:	2 - 100 µg/L (ppb) Cl ₂ free
Number of tests:	100
Parameter:	Free Chlorine
What's included?:	Includes 1 bottle Free Chlorine Reagent 1, 1 bottle Free Chlorine Reagent 2, 2x 16 mm Sample Vials, Vial Stand

What's included?

Includes 1 bottle Free Chlorine Reagent 1, 1 bottle Free Chlorine Reagent 2, 2x 16 mm Sample Vials, Vial Stand