



ULR Total Chlorine Fluorescence Reagent Refill Kit, 3-100 µg/L (ppb), 100 Tests

Product #:

34252001

Hazardous

AED Price:

Contact Hach



Refill Reagents for ULR Total Chlorine Test Kit 34252000.

Specifications

Measuring range:	3 - 100 µg/L (ppb) Cl ₂ total
Number of tests:	100
Parameter:	Total Chlorine
What's included?:	Includes 1 bottle Total Chlorine Reagent 1, 1 bottle Total Chlorine Reagent 2

What's included?

Includes 1 bottle Total Chlorine Reagent 1, 1 bottle Total Chlorine Reagent 2