



ULR Sulfite Fluorescence Reagent Refill Kit, 6-500 µg/L (ppb), 100 Tests

Product #: 34250001
AED Price: Contact Hach

Refill Reagents for ULR Sulfite Test Kit 34250000.

Specifications

Measuring range:	6 - 500 µg/L (ppb) SO ₃
Number of tests:	100
Parameter:	Sulfite
What's included?:	Includes 1 bottle Sulfite Reagent 1, 1 bottle Sulfite Reagent 2

What's included?

Includes 1 bottle Sulfite Reagent 1, 1 bottle Sulfite Reagent 2