



## EDTA solution for EZ4043 Total Hardness Analyser, 0.005 N, 2 L

Product #:	APPC4043-04
AED Price:	Contact Hach

### Ready-to-use Reagents for EZ Series Online Analysers

#### Ready-to-use

Hach EZ reagents are ready-to use with EZ Series Online Analysers. Since there is no additional preparation needed this eliminates an additional error source when using the EZ instruments and reagents system.

#### Confidence in your results

Reduce the sources of errors in your workflow to increase the confidence in your results.

#### Less complexity

Using ready-to-use EZ reagent kits from Hach can reduce the number of handling steps around your analytical equipment significantly which allows less hands-on-time and complexity.

---

### Specifications

Concentration:	0.005 N
Parameter:	Total Hardness
Quantity:	2 L