



Reagent Kit for EZ4043 Total Hardness Analyser (acid, buffer, colour and EDTA solution), 0.25-10 mg/L CaCO₃

Product #:

APPC4043KTO-EU **Hazardous**

AED Price:

Contact Hach



Ready-to-use Reagents for EZ Series Online Analysers

Ready-to-use

Hach EZ reagents are ready-to use with EZ Series Online Analysers. Since there is no additional preparation needed this eliminates an additional error source when using the EZ instruments and reagents system.

Confidence in your results

Reduce the sources of errors in your workflow to increase the confidence in your results.

Less complexity

Using ready-to-use EZ reagent kits from Hach can reduce the number of handling steps around your analytical equipment significantly which allows less hands-on-time and complexity.

Specifications

Parameter:	Total Hardness
What's included?:	APPC4041-01-EU Acid Solution
	APPC4041-02-EU Buffer Solution
	APPC4041-03 Colour Solution
	APPC4043-04 EDTA Solution, 0.005 N

What's included?

APPC4041-01-EU Acid Solution APPC4041-02-EU Buffer Solution APPC4041-03 Colour Solution APPC4043-04 EDTA Solution, 0.005 N